

Walkup MS Report



Data File	SE-455 A1(AC)_Neg_5mins_MS_07310.d	Sample Name	SE-455 A1(AC)
Sample Type	Sample	Position	P1-A4
Instrument Name	6545 QTof	User Name	Susana Estopina-Duran
Acq Method	Neg_5mins_MS.m	Acquired Time	8/29/2019 9:06:22 AM
IRM Calibration Status	Success	DA Method	Neg_5mins_MS.m
Comment			

Sample Group		Info.	
Walkup Sample Description		Walkup Method	Neg_5Mins_C18
Formula	C19H18O2	Walkup Method Description	Negative mode ionization using C18 column chromatography
Stream Name	LC 1	Acquisition SW Version	6200 series TOF/6500 series Q-TOF B.09.00 (B9044.0)




Figure 1: Base peak chromatogram

User Chromatogram Peak List

RT (min)	Area	Area %	Area Sum (%)	Base Peak (m/z)	Width (min)
4.01	4974621	11.38	3.05	1033.9836	0.140
4.22	5046342	11.54	3.10	1033.9839	0.120
4.40	3481346	7.96	2.14	1033.9846	0.110
4.63	32682044	74.76	20.05	277.1231	0.140
4.97	43713165	100.00	26.82	265.1567	0.180
5.10	18186665	41.60	11.16	311.2007	0.090
5.21	25832215	59.09	15.85	305.2146	0.170
5.51	6965695	15.94	4.27	1033.9854	0.090
5.66	21535523	49.27	13.21	265.1481	0.190
6.08	547375	1.25	0.34	1033.9862	0.160

Compound Table

Compound Label	RT (min)	Observed mass (m/z)	Neutral observed mass (Da)	Theoretical mass (Da)	Mass error (ppm)	Isotope match score (%)
Cpd 1: C19 H18 O2	6.34	277.1229	278.1302	278.1307	-1.82	98.96

Mass errors of between -5.00 and 5.00 ppm with isotope match scores above 60% are considered confirmation of molecular formulae

Walkup MS Report



Compound specific information




Figure: Extracted ion chromatogram (EIC) of compound.




Figure: Full range view of Compound spectra and potential adducts.




Figure: Zoomed Compound spectra view
(red boxes indicating expected theoretical isotope spacing and abundance)

Compound isotope peak List

m/z	z	Abund	Formula	Ion
277.1229	1	7815.7	C ₁₉ H ₁₈ O ₂	(M-H)-
278.1263	1	1615.3	C ₁₉ H ₁₈ O ₂	(M-H)-
279.1289	1	222.4	C ₁₉ H ₁₈ O ₂	(M-H)-

--- End Of Report ---