

Single Mass Analysis

Tolerance = 5.0 PPM / DBE: min = -1.5, max = 100.0
Element prediction: Off
Number of isotope peaks used for i-FIT = 3

Monoisotopic Mass, Odd and Even Electron Ions
2 formula(e) evaluated with 1 results within limits (up to 50 closest results for each mass)
Elements Used:
C: 0-17 H: 0-16 O: 0-1 F: 0-1

Minimum:				-1.5				
Maximum:	5.0	5.0		100.0				
Mass	Calc. Mass	mDa	PPM	DBE	i-FIT	Norm	Conf (%)	Formula
255.1187	255.1185	0.2	0.8	9.5	1153.9	n/a	n/a	C17 H16 O F