

Single Mass Analysis

Tolerance = 5.0 PPM / DBE: min = -1.5, max = 100.0

Element prediction: Off

Number of isotope peaks used for i-FIT = 3

Monoisotopic Mass, Odd and Even Electron Ions
63 formula(e) evaluated with 1 results within limits (up to 50 closest results for each mass)
Elements Used:
C: 0-16 1H: 0-13 O: 0-2 F: 0-1 2H: 0-1

Minimum:				-1.5					
Maximum:	5.0	5.0		100.0					
Mass	Calc. Mass	mDa	PPM	DBE	i-FIT	Norm	Conf (%)	Formula	
258.1043	258.1041	0.2	0.8	9.5	1156.0	n/a	n/a	C16 1H13 O2 F 2H	