

Single Mass Analysis

Tolerance = 5.0 PPM / DBE: min = -1.5, max = 100.0
Element prediction: Off
Number of isotope peaks used for i-FIT = 3

Monoisotopic Mass, Odd and Even Electron Ions
7 formula(e) evaluated with 1 results within limits (up to 50 closest results for each mass)
Elements Used:
C: 0-16 H: 0-13 O: 0-2 F: 0-2

Minimum:				-1.5				
Maximum:	5.0	5.0		100.0				
Mass	Calc. Mass	mDa	PPM	DBE	i-FIT	Norm	Conf (%)	Formula
275.0897	275.0884	1.3	4.7	9.5	1342.8	n/a	n/a	C16 H13 O2 F2