

Single Mass Analysis

Tolerance = 5.0 PPM / DBE: min = -1.5, max = 100.0
Element prediction: Off
Number of isotope peaks used for i-FIT = 3

Monoisotopic Mass, Odd and Even Electron Ions
2 formula(e) evaluated with 1 results within limits (up to 50 closest results for each mass)
Elements Used:
C: 0-23 H: 0-19 N: 0-2 O: 0-1

Minimum:				-1.5				
Maximum:	5.0	5.0		100.0				
Mass	Calc. Mass	mDa	PPM	DBE	i-FIT	Norm	Conf (%)	Formula
339.1498	339.1497	0.1	0.3	15.5	1245.6	n/a	n/a	C23 H19 N2 O