

Single Mass Analysis
Tolerance = 10.0 PPM / DBE: min = -1.5, max = 100.0
Element prediction: Off
Number of isotope peaks used for i-FIT = 3

Monoisotopic Mass, Odd and Even Electron Ions
3 formula(e) evaluated with 1 results within limits (up to 50 closest results for each mass)
Elements Used:
C: 0-16 H: 0-14 N: 0-1 O: 0-2

Minimum:				-1.5				
Maximum:		5.0	10.0	100.0				
Mass	Calc. Mass	mDa	PPM	DBE	i-FIT	Norm	Conf (%)	Formula
252.1033	252.1025	0.8	3.2	10.5	1311.0	n/a	n/a	C16 H14 N O2