

Single Mass Analysis

Tolerance = 5.0 PPM / DBE: min = -1.5, max = 100.0

Element prediction: Off

Number of isotope peaks used for i-FIT = 3

Monoisotopic Mass, Odd and Even Electron Ions
3 formula(e) evaluated with 1 results within limits (up to 50 closest results for each mass)
Elements Used:
C: 0-15 H: 0-15 O: 0-3

Minimum:				-1.5				
Maximum:	5.0	5.0		100.0				
Mass	Calc. Mass	mDa	PPM	DBE	i-FIT	Norm	Conf (%)	Formula
243.1020	243.1021	-0.1	-0.4	8.5	1129.5	n/a	n/a	C15 H15 O3