

Single Mass Analysis

Tolerance = 5.0 PPM / DBE: min = -1.5, max = 100.0

Element prediction: Off

Number of isotope peaks used for i-FIT = 3

Monoisotopic Mass, Odd and Even Electron Ions
8 formula(e) evaluated with 1 results within limits (up to 50 closest results for each mass)
Elements Used:
C: 0-18 H: 0-16 O: 0-4 F: 0-1

Minimum:				-1.5				
Maximum:	5.0	5.0		100.0				
Mass	Calc. Mass	mDa	PPM	DBE	i-FIT	Norm	Conf (%)	Formula
315.1031	315.1033	-0.2	-0.6	10.5	895.1	n/a	n/a	C18 H16 O4 F