

Pure zwitterionic surfactants in choline chloride:glycerol

Sample	Surfactant	Solvent	Conc / mM
NC201	h-C12-PC	1:2 d-choline chloride:d-glycerol	31.5±1.2
NC202	h-C12-PC	1:2 d-choline chloride:d-glycerol	68.1±2.4
NC203	h-C12-PC	1:2 d-choline chloride:d-glycerol	134±3
NC204	h-C12-PC	1:2 d-choline chloride:d-glycerol	308±2
NC205	h-C12-PC	1:2 h/d-choline chloride:h/d-glycerol	31.5±1.2
NC206	h-C12-PC	1:2 h/d-choline chloride:h/d-glycerol	68.1±2.4
NC207	h-C12-PC	1:2 h/d-choline chloride:h/d-glycerol	134±3
NC208	h-C12-PC	1:2 h/d-choline chloride:h/d-glycerol	308±2
NC209	h-C12-PC	1:2 h-choline chloride:d-glycerol	31.5±1.2
NC210	h-C12-PC	1:2 h-choline chloride:d-glycerol	68.1±2.4
NC211	h-C12-PC	1:2 h-choline chloride:d-glycerol	134±3
NC212	h-C12-PC	1:2 h-choline chloride:d-glycerol	308±2
NS201	h-SB3-12	1:2 d-choline chloride:d-glycerol	28.8±0.9
NS202	h-SB3-12	1:2 d-choline chloride:d-glycerol	64.3±1.8
NS203	h-SB3-12	1:2 d-choline chloride:d-glycerol	151±5
NS204	h-SB3-12	1:2 d-choline chloride:d-glycerol	312±6
NS205	h-SB3-12	1:2 h/d-choline chloride:h/d-glycerol	28.8±0.9
NS206	h-SB3-12	1:2 h/d-choline chloride:h/d-glycerol	64.3±1.8
NS207	h-SB3-12	1:2 h/d-choline chloride:h/d-glycerol	151±5
NS208	h-SB3-12	1:2 h/d-choline chloride:h/d-glycerol	312±6
NS209	h-SB3-12	1:2 h-choline chloride:d-glycerol	28.8±0.9
NS210	h-SB3-12	1:2 h-choline chloride:d-glycerol	64.3±1.8
NS211	h-SB3-12	1:2 h-choline chloride:d-glycerol	151±5
NS212	h-SB3-12	1:2 h-choline chloride:d-glycerol	312±6
NS401	h-SB3-14	1:2 d-choline chloride:d-glycerol	13.6±1.1
NS402	h-SB3-14	1:2 d-choline chloride:d-glycerol	30.0±3.1
NS403	h-SB3-14	1:2 d-choline chloride:d-glycerol	63.1±1.6
NS404	h-SB3-14	1:2 d-choline chloride:d-glycerol	166±4
NS405	h-SB3-14	1:2 h/d-choline chloride:h/d-glycerol	13.6±1.1
NS406	h-SB3-14	1:2 h/d-choline chloride:h/d-glycerol	30.0±3.1
NS407	h-SB3-14	1:2 h/d-choline chloride:h/d-glycerol	63.1±1.6
NS408	h-SB3-14	1:2 h/d-choline chloride:h/d-glycerol	166±4
NS409	h-SB3-14	1:2 h-choline chloride:d-glycerol	13.6±1.1
NS410	h-SB3-14	1:2 h-choline chloride:d-glycerol	30.0±3.1
NS411	h-SB3-14	1:2 h-choline chloride:d-glycerol	63.1±1.6
NS412	h-SB3-14	1:2 h-choline chloride:d-glycerol	166±4

Mixtures of zwitterionic surfactants in choline chloride:glycerol

Sample	Ratio h-SB3-12:h-C12-PC	Solvent	Conc / mM
Mx01	0.2:0.8	1:2 d-choline chloride:d-glycerol	72.3±1.2
Mx02	0.2:0.8	1:2 d-choline chloride:d-glycerol	185±3
Mx03	0.35:0.65	1:2 d-choline chloride:d-glycerol	72.3±1.2
Mx04	0.35:0.65	1:2 d-choline chloride:d-glycerol	185±3
Mx05	0.5:0.5	1:2 d-choline chloride:d-glycerol	72.3±1.2
Mx06	0.5:0.5	1:2 d-choline chloride:d-glycerol	185±3
Mx07	0.65:0.35	1:2 d-choline chloride:d-glycerol	72.3±1.2
Mx08	0.65:0.35	1:2 d-choline chloride:d-glycerol	185±3
Mx09	0.8:0.2	1:2 d-choline chloride:d-glycerol	72.3±1.2
Mx10	0.8:0.2	1:2 d-choline chloride:d-glycerol	185±3